

EN

GD4210

MAVIG



GD4210 Rear view

Article-No.:	GD4210	Monitor support for one display up to 19", with VESA 75/100 adapter
	GD4210-L	Monitor support for one display up to 21/24", with VESA 75/100 adapter
	GD4210-XL	Monitor support for one display up to 27" (some up to 31"), with VESA 75/100 adapter


Made in Germany

GD4210 – Single Monitor Support* for Portegra2 Systems

The GD4210 can easily hold one monitor up to 27" (some up to 31"). The large curved handlebar protects the monitor and provides a means for quick & easy positioning. Combined with the mobility and safety of the Portegra2 suspension arm, the GD4210 offers the user a worry-free and flexible application.

- **Application** Monitor support with handle, VESA 75/100 adapter for one display/monitor, external cabling
- **Load / Weight**

GD4210	max. 15.5 kg (33 lbs) / 2.5 kg
GD4210-L	max. 15.4 kg (33 lbs) / 2.6 kg
GD4210-XL	max. 15.2 kg (33 lbs) / 2.8 kg
- **VESA adapter** VESA 75/100 compatible
- **Colouring** RAL9002
- **Mounting** Portegra2 systems for Stationary / ceiling track mounting
- **System range** Extension arm 750 mm / 950 mm
Spring arm 910 mm



* (Dimensions for L and XL versions)

* Monitor support without extension / spring arm, cabling and monitor.

MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДІКА“

**Уповноважений
представник в Україні**
01014, Україна, м. Київ, вул.
Звіринська, 63
+38 (044) 228 84 17
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua