

EN

OT50001I

MAVIG



STEA-OT5 (accessory) Disposable, sterile cover for the whole shield

Article Codes: OT50001I Shield complete with Portegra2 Suspension Arm (extension arm length 75 cm)
E-OT25B05I Replacement Shield and Connection Bar

IEC 60601-1

Made in Germany




OT50001I – X-Ray Protective Shield

MAVIG's products are designed to offer solid protection and optimized functionality. The high quality acrylic shield allows for an excellent transmittance and unimpeded view. All edges and corners have been rounded and shatter-proof materials are used to minimize risk of injury.

MAVIG's patented Portegra2 suspension arm provides its industry known safety features and movement flexibility. The centrally guided shield with a special patient cut-out can be easily positioned over the patient and quickly removed in case of an emergency.

- **Lead Acrylic Shield** WxH: 60 x 76 cm (23.6 x 29.9 in)
- **Lead Equivalent** Pb 0.50 mm
- **OT50001I** Extension / Spring Arm
Lengths: 75 cm / 91 cm (29.5 in / 35.8 in)
Colour: RAL 9002



MAVIG GmbH

X-Ray Protection and Medical Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com
ТОВ „ОНКОМЕДИКА“
Уповноважений представник в Україні
01014, Україна, м. Київ, вул. Звіринецька, 63
+38 (044) 228 84 17
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua