


# ST – X-Ray Protective Drapes

## ST-FS5AMM



## ST-RZ5AMM



|                |                  |   |
|----------------|------------------|---|
| <b>Models:</b> | <b>ST-FS5AMM</b> | X-ray protective drape with cutout for femoral access * |
|                | STEA-FSAM        | 50x sterile covers for X-ray protective drape ST-FS5AMM |
|                | <b>ST-RZ5AMM</b> | X-ray protective drape for radial access *              |
|                | STEA-RZAM        | 50x sterile covers for X-ray protective drape ST-RZ5AMM |



Made in Germany

DIN EN 61331-1  
IEC 61331-1

## X-Ray Protective Drapes: ST-FS5AMM & ST-RZ5AMM

Two uniquely formed, protective drapes have been designed to optimally protect against radiation emitting from the patient's body during both femoral and lateral access procedures.

Our covers have been ergonomically formed to provide excellent protection, while being as compact as possible for patient comfort. Special attention has also been made to their application, such as the special cutout of the femoral cover for the femoral port, and the gentle slope of the right side of both covers for unobstructed imaging.

Combined with our dedicated sterile covers, the protective drapes can be placed on top of the sterile patient drape. This allows for quick placement or removal, as well as constant reference to the cover's position in regards to the image field. The sterile covers also provide an economic and environmental advantage to common disposable drapes.

### ST-FS5AMM

- Drape with cutout for femoral access
- Lead equivalent: Pb 0.50 mm
- Size: 75 x 36 cm / 29.5 x 14.2 in
- Weight: 1.31 kg / 2.89 lbs

### ST-RZ5AMM

- Drape for radial access
- Lead equivalent: Pb 0.50 mm
- Size: 75 x 40 cm / 29.5 x 15.7 in
- Weight: 1.56 kg / 3.44 lbs

\* including one box of 50 sterile covers

### MAVIG GmbH

X-Ray Protection and Medical  
Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua