


PORTEGRA<sup>2</sup>

## TriPin Ceiling Column

## Stationary Mounting

Two Electrical Connections  
with 3-pole sliding contacts



Installation Example

<b>Models:</b>	<b>TS2320</b>	TriPin column for stationary mounting, length 1000 mm, with safety steel rope
	<b>TS2321</b>	TriPin column for stationary mounting, length 850 mm, with safety steel rope
	<b>TS2322</b>	TriPin column for stationary mounting, length 580 mm, with safety steel rope

Made in Germany

IEC 60601-1

## TS232x – Portegra2 TriPin Ceiling Column

The innovative TriPin ceiling column series expands on Portegra2's successful slim design and flexibility to provide a cost-effective way to manage three devices within a compact space. Portegra2's signature 360° mounting pin still guarantees the highest degree of moveability and positioning.

MAVIG TriPin columns are equipped standard with two electrical connections for fixed mounted systems.

- **Three mounting posts** (two upper, one lower)
- **Both upper posts offer a 235° rotation; Lower pin 360°**
- **Two electrical connection with 3-pole sliding contacts** ⚡
- **Maximum 18.0 kg (39.7 lbs) net load for each mounting post (depending on length)**
- **Different lengths** (580 mm, 850 mm, or 1000 mm) **to match room height**
- **Compatible with existing Portegra2 mounting plates and suspension arms**

Special lengths upon request: TriPin column stationary: minimal length = 530 mm

Ceiling column lengths are measured from the top of the mounting flange to the bottom of the column structure.

**MAVIG GmbH**

X-Ray Protection and Medical  
Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com

**ТОВ „ОНКОМЕДИКА“**

**Уповноважений  
представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua