

EN

BR119

Designs Silver
Black

MAVIG



Silver



Black



Model: BR119 Glasses with front X-ray protection of 0.75 mmPb, universal size, weight approx. 65 g

CE 2895

BR119 – X-Ray Protective Glasses

The X-ray protective glasses BR119 have curved glasses and a front protection of 0.75 mmPb. The frame is made of an extremely durable plastic that ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR119 comes with a storage case and a soft google strap.

MAVIG uses only protective lenses of the highest quality. This guaranties the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

The quality features of our glasses make them perfect for the use in the medical environment:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- For this model prescription lenses are available: single, bifocal, and progressive
- Conformity: DIN EN 61331-1:2016 / DIN EN 61331-3:2016

MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems

Stahlgruberring 5
81829 Munich

Germany
PO Box 82 03 62

81803 Munich

Germany

Phone +49 (0) 89 / 420 96 - 0

Fax: + 49 (0) 89 / 420 96 - 200

e-Mail info@mavig.com

www.mavig.com

ТОВ „ОНКОМЕДИКА“

Уповноважений

представник в Україні

01014, Україна, м. Київ, вул.

Звіринецька, 63

+38 (044) 228 84 17

+38 (097) 796 71 30

+38 (097) 959 10 89

info@oncomedita.ua

www.oncomedita.ua

www.mavig.com.ua