

EN

# BRV500 BRV501



Model BRV500



Model BRV501

**Models:** **BRV500** Continuous protection from forehead to the chin  
**BRV501** Protection down to the cheek bones with a cut-out for the nose

## BRV500 & BRV501 – X-Ray Protective Visors

Especially designed for the protection of medical professionals during radiological procedures. Protection of the eye lens against scattered X-rays. The visors provide a high reduction in radiation exposure to the eyes and most of the facial area (see table).

**The quality features of our visors make them ideal for the use in the medical environment**

- Pb Value (Weight): BRV500 0.10 mmPb (approx. 520 g)  
BRV501 0.10 mmPb (approx. 480 g)
- Universal size / adjustable to fit different head sizes
- Multiple adjustments to conform to individual facial features
- Easy to set up for fitting the head of the wearer by simply pressing and turning the pivoting mechanism to prevent slipping
- Comfortable, exchangeable, and washable pads
- Curved lead acrylic glass with 0.10 mmPb
- Low allergy risk due to carefully chosen materials

X-Ray Tube Voltage	Attenuation Scattered Radiation
50 kV	96.9 %
60 kV	94.4 %
80 kV	86.3 %
100 kV	80.0 %
120 kV	75.0 %

**MAVIG GmbH**  
X-Ray Protection and Medical Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com  
**ТОВ „ОНКОМЕДІКА”**  
**Уповноважений представник в Україні**  
01014, Україна, м. Київ, вул. Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 796 71 30  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua