

# EN OR HD-Camera System

MAVIG



**Camera Housing**  
integrated in the OR-lamp

**Camera - Remote Control**  
included



Combined  
with LED6 MC Lamp

## Features

- 30-fold optical zoom, 12-fold digital zoom
- automatic/manual focus
- Image rotation without end stops
- Automatic white balance
- HD-resolution: 1920 x 1080 pixels
- automatic/manual brightness control
- Freeze frame function
- Remote controlled

**Model:** KAM02

Wireless HD-camera for usage with LED6 MC lamp (item LE7069111KVS), control unit and remote control included

Made in Germany



## Wireless HD-Camera System for OR Lamps

The camera system with a live video image in HD resolution and high light sensitivity (minimum illumination 1.4 lx). The brilliant colour reproduction, even in difficult lighting conditions, is especially suitable for recording operations in the operating theatre for documentation and for training/education.

Thanks to the lamp design, the camera is always precisely aligned at the centre of the illuminated field. Power is supplied directly via the lamp. The camera also has an electronic device to control focus, zoom and lens iris. All available functions can also be operated via the control unit.

- **Camera lens system** 30-fold optical zoom, 12-fold digital zoom  
F = 4.3 to 129 mm (F 1.6 - 4.7, auto-focus)  
Minimum 70 mm (wide) to 140 mm (tele)  
1920 x 1080 / approx. 2.000.000 pixels, aspect ratio 16:9  
1080p – 50 Hz, 60 Hz / 1/2,8-type CMOS  
50 dB
- **Object distance** 1x HDMI (incl. 10 m HDMI cable)
- **Resolution / Image points** 80 x 190 mm (Ø, length) / 900 g (HD camera wireless)  
in acc. with FCC class A
- **Video signal / Image sensor** in accordance with MDR Regulation (EU) 2017/745
- **S/N ratio** CE
- **Video outputs**
- **Dimensions / Weight**
- **Interference radiation**
- **Protection and approval**
- **Medical Product - (EU) 2017/745, Usage in accordance with DIN VDE 0100-710**

## MAVIG GmbH

X-Ray Protection and Medical Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63

+38 (044) 228 84 17  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua