

EN

RA611

with breathable insert



Included Outlast
Temperature regulation

MAVIG

RA612

closed cap



Included Outlast
Temperature regulation



Models: RA611
RA612

Radiation protection cap, with breathable insert ensuring very good heat dissipation
Radiation protection cap, closed cap shielding the cranium frontally, on the sides and completely from above (limits heat dissipation)

Made in Germany

Radiation Protection Head Wear – RA611 & RA612

Working in interventional radiology can entail risks. Currently, the likelihood of medical personnel being exposed to radiation and developing a brain tumour is the subject of lively debate among experts.

Without effective on-site upper table scatter radiation protection, the cranium may well be exposed to high levels of radiation, depending on the procedure and fluoroscopy time.

MAVIG, therefore, has developed two different radiation protective head coverings to meet the individual needs of the user. The respective models differ in the matching of the desired protection area and heat elimination.

- **Effective protection for the cranium** (optional: 0.25 mm Pb, 0.35 mm Pb, or 0.50 mm Pb)
- **Outer material: ComforTex® HPMF**
- **Corresponds to work-related physiological aspects**
- **Perfect balance of radiation protection, weight and comfort**
- **NEW: Always with Outlast® material**
- **Design-protected, universal-size cuts**
- **Optimal adjustment to the respective head circumference and head shape**
- **Sizes: suitable for head circumferences from 54 - 63 cm**

MAVIG GmbH
X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДИКА“
Уповноважений
представник в Україні
01014, Україна, м. Київ, вул.
Звіринецька, 63

+38 (044) 228 84 17
+38 (097) 796 71 30
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua