

Basic RA640

FRONT PROTECTION APRON

The basic RA640 model has proven effective for years and is used as front protection usually for short and seldom wearing times in conventional X-ray situations without fluoroscopy.

Protection/Lead Equivalent Standard Collection

NovaLite, Leadfree or Standard lead



Front:

0.25 mm Pb, 0.35 mm Pb or 0.50 mm Pb

Back generally:
0.25 mm Pb

- ✓ Ideal fit, maximum comfort for the wearer
- ✓ Ideal for short and infrequent wear
- ✓ Soft shoulder pads as standard equipment
- ✓ Various closure types to choose from
- ✓ Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Sizes*, recommended for:

XS	Chest girth ≤ 84 cm	Waist/hip girth ≤ 91 cm
S	Chest girth approx. 85 - 95 cm	Waist/hip girth approx. 92 - 102 cm
M	Chest girth approx. 96 - 106 cm	Waist/hip girth approx. 103 - 112 cm
L	Chest girth approx. 107 - 117 cm	Waist/hip girth approx. 113 - 123 cm
XL	Chest girth approx. 118 - 128 cm	Waist/hip girth approx. 124 - 134 cm
XXL	Chest girth approx. 129 - 139 cm	Waist/hip girth approx. 135 - 145 cm
3XL	Chest girth approx. 140 - 150 cm	Waist/hip girth approx. 146 - 156 cm

* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.







Outlast Temperature Regulation

Different Closure Types – You can choose ...


We offer different closure systems for this series of aprons. Choose your favourite:




Snap fastener closure
Elastic band with snap fasteners for closing in the front.



Crossband closure
Crossbands with hook-and-loop fastening for closing in the back.



Strap closure
Straps with snap buckle closure for closing in the front.




Snap buckle closure
Two webbing straps with snap buckles for closing in the back.



Back wings

For model RA640 it is also possible to choose whether a protection of 0.25 mm Pb should be maintained in the back panels or not, in order to reduce the weight of the apron.



Outer Material

Colours
(see front page)



Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid

Sizes

from XS to 3XL
(see front page)



Length / Body height:
90 cm / 150 - 160 cm
100 cm / 160 - 170 cm
110 cm / 170 - 180 cm
120 cm / 175 - 190 cm
130 cm / 185 - 195 cm
140 cm / 195 - 205 cm

Optional



On request:
Outlast
Sleeves
Breast pocket
Embroidery
Made-to-measure protection

Accessories/Replacement

Bolero
Thyroid/sternum protection



Crossband closure
(Adjustment via
snap fastener)