

# Comfort RA660

## FRONT PROTECTION APRON

Our RA660 range offers up-to-date radiation protection clothing which accompanies your active working day and provides outstanding service. The extraordinarily good wearing comfort is supported by the wide stretch insert integrated in the back panels.

### Protection/Lead Equivalent Standard Collection

NovaLite, Leadfree  
or Standard lead



**Front:**  
0.25 mm Pb, 0.35 mm Pb  
or 0.50 mm Pb

**Back generally:**  
0.25 mm Pb

### Protection/Lead Equivalent Levio Edition

NovaLite, Leadfree  
or Standard lead



**Total value Levio-Panel:**  
0.42 mm Pb or 0.50 mm Pb  
(outside 0.25 mm Pb)

**Back generally:**  
0.25 mm Pb

- ✔ Suitable for a wide range of applications with only frontal radiation exposure
- ✔ Ideal fit, maximum comfort for the wearer
- ✔ Suitable for long wearing times
- ✔ Optionally "without protection" in the back, for a reduced weight of the apron
- ✔ Standard equipment includes shoulder pads
- ✔ Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Sizes\*, recommended for:

XS	Chest girth ≤ 84 cm	Waist/hip girth ≤ 91 cm
S	Chest girth approx. 85 - 95 cm	Waist/hip girth approx. 92 - 102 cm
M	Chest girth approx. 96 - 106 cm	Waist/hip girth approx. 103 - 112 cm
L	Chest girth approx. 107 - 117 cm	Waist/hip girth approx. 113 - 123 cm
XL	Chest girth approx. 118 - 128 cm	Waist/hip girth approx. 124 - 134 cm
XXL	Chest girth approx. 129 - 139 cm	Waist/hip girth approx. 135 - 145 cm
3XL	Chest girth approx. 140 - 150 cm	Waist/hip girth approx. 146 - 156 cm

\* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.



## Sales Champion

Our customers report that this model is their favourite front apron. This is because it can be put on and taken off easily, fits the body well, stays in place well and provides excellent shoulder relief.


## Closure System

MAVIG uses only top quality, low-wear hook-and-loop fasteners.

Levio®   
Optimised Protection Concept



Outlast Temperature Regulation



### Stretch insert

The stretch insert running across the entire width of the back ensures a good fit and shoulder relief even when the X-ray protective clothing is worn for long periods.

### Higher Protection – Levio Edition

The area with the reinforced protective material is clearly recognisable by its colour.

### User Friendly

The sophisticated design of the RA660 front protection apron makes it easy and quick to put on or take off.

### Special Feature

The stretch insert offers ideal relief especially for long periods of wear.



### Outer Material

Colours  
(see front page)



Optional:  
ComforTex® HPMF,  
ComforTex® HPMF Hybrid

### Sizes

from XS to 3XL  
(see front page)



Lengths / Body height:  
90 cm / 150 - 160 cm  
100 cm / 160 - 170 cm  
110 cm / 170 - 180 cm  
120 cm / 175 - 190 cm  
130 cm / 185 - 195 cm  
140 cm / 195 - 205 cm

### Optional



On request:  
Outlast  
Sleeves  
Breast pocket  
Embroidery  
Made-to-measure protection

### Accessories

Bolero  
Thyroid/sternum protection

