


EN  
RD614

MAVIG




**Note:** In order to meet hygiene requirements, we recommend that you use our washable hygienic cover made from special ComforTex® MF material "Titan" colour.



With interchangeable hook-and-loop fastener



Children's size in Happy design



Available optional hygienic cover

<b>Models:</b>	<b>RD6145E</b>	Adult patient thyroid and sternum protection
	<b>RD6145K</b>	Children patient thyroid and sternum protection

Made in Germany

CE 0302

## RD614 – Thyroid and Sternum Protection for Patients

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection. The MAVIG thyroid sternum protection for patients is optimally cut and has the best fit, and, of course, it also takes all radiation protection aspects into account.

An uncomplicated magnetic or hook-and-loop fastener, which can also be replaced as an option, provides sufficient scope for adjustment.

- **Protective Material** Standard lead or Leadfree
- **Lead Equivalent** 0.50 mm Pb
- **Magnetic/hook-and-loop fastener**

Adults	Standard size fits neck sizes of 32 - 45 cm
Children	Standard size fits neck sizes of 28 - 36 cm
- **Interchangeable hook-and-loop fastener**

Adults	Standard size fits neck sizes of 37 - 50 cm
Children	Standard size fits neck sizes of 28 - 36 cm
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)*

### MAVIG GmbH

X-Ray Protection and Medical Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 796 71 30  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua