

EN

RD642

MAVIG



**Models:** **RD6425E** Adult patient X-ray apron for X-rays in a seated position  
**RD6425K** Children patient X-ray apron for X-rays in a seated position

Made in Germany

CE 0302

## RD642 – Patient X-Ray Apron for X-Rays in a Seated Position

The MAVIG patient protective apron RD642, which is used for dental X-rays in a seated position, is flexible and comfortable to wear.

The high-quality hook-and-loop fasteners are conveniently placed in the shoulder area, making it precise, secure and quick to put on and remove.

- **Protective Material** Standard lead or Leadfree
- **Lead Equivalent** 0.50 mm Pb (front)
- **Size**

Adult	Front: 80.0 x 60.0 cm (H x W)
Children	Front: 60.0 x 45.0 cm (H x W)
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy (ComforTex® HPMF only)*

**MAVIG GmbH**  
 X-Ray Protection and  
 Medical Suspension Systems  
 Stahlgruberring 5  
 81829 Munich  
 Germany  
 PO Box 82 03 62  
 81803 Munich  
 Germany  
 Phone +49 (0) 89 / 420 96 - 0  
 Fax: +49 (0) 89 / 420 96 - 200  
 e-Mail info@mavig.com  
 www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
 01014, Україна, м. Київ, вул.  
 Звіринецька, 63  
 + 38 (044) 228 84 17  
 + 38 (097) 796 71 30  
 + 38 (097) 959 10 89  
 info@oncomedica.ua  
 www.oncomedica.ua  
 www.mavig.com.ua