

EN RD644

MAVIG



CAN BE USED ON BOTH SIDES

The model RD644 is designed as a reversible X-ray apron for patients.



Models:
RD644SE
RD644SK

Adult patient X-ray apron with thyroid protection, for X-rays in a seated position
Children patient X-ray apron with thyroid protection, for X-rays in a seated position

Made in Germany

CE 0302

RD644 – Patient X-Ray Apron with Thyroid Protection

Just like the classic dental apron (RD642), this model offers high flexibility and wearing comfort. In addition to the large-area protection up to the gonads, the thyroid gland is also covered when taking X-rays of the seated patient.

The practical hook-and-loop fastener at the back of the neck on the thyroid protection collar ensures a snug, secure and quick fit. Loops on the back allow easy attachment to the dentist's chair.

- **Protective Material** Standard lead or Leadfree
- **Lead Equivalent** 0.50 mm Pb (front)
- **Size**
 - Adult Front: 80.0 x 60.0 cm (H x W)
 - Children Front: 60.0 x 45.0 cm (H x W)
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Hybrid, available in the colours:
 - Titan, Blue Admiral, Regatta, Curacao,
 - Lagoon, Mandarin, Orchidee, Baccara,
 - or design: Happy (ComforTex® HPMF only)
 - and two-colour combinations:
 - Curacao with Orchidee or
 - Blue Admiral with Baccara

MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5 81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДІКА“ Уповноважений представник в Україні

01014, Україна, м. Київ, вул.
Звіринська, 63
+ 38 (044) 228 84 17
+ 38 (097) 796 71 30
+ 38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua