



Model: RP2705 Breast protection for female patients during mammographic radiation procedures

Made in Germany


CE 0302

RP270 – Breast Protection for Women

The mammary gland shield efficiently protects your patients during mammograms.

The shield consists of two equal parts, each 18 x 40 cm (W x H), connected by a lengthadjustable hook-and-loop fastener for comfortable donning and easy adjustment.

- **Protective Material** Novalite radiation protective material
- **Lead Equivalent** 0.50 mm Pb
- **Size** Universal
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or designs: Happy
- **Edging** All edges are covered with a seem band.



MAVIG GmbH

X-Ray Protection and Medical Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДІКА“
Уповноважений представник в Україні
01014, Україна, м. Київ, вул. Звіринецька, 63
+38 (044) 228 84 17
+38 (097) 796 71 30
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua