

EN
RP271

MAVIG



Size Adjustment From Behind

Model: RP27156231 Mammary Gland Cover for Children

Made in Germany


CE 0302

RP271 – X-Ray Protection for Paediatrics

The patient protection product for paediatrics reliably protects the mammary glands of your little patients during X-rays in medical diagnostics. The children's chest protector can be conveniently adjusted to the respective size required.

- ✓ Reduction of radiation exposure in the area of the mammary glands by about 80%
- ✓ Hygienically clean, skin-friendly outer cover
- ✓ Reusable, not a disposable product
- ✓ Comfortable, non-pressing neoprene stoppers
- ✓ Practical adjustment

- **Protective Material** Novalite radiation protective material
- **Lead Equivalent** 0.50 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF Plus available in design Happy
- **Size** Universal, easily adjustable



MAVIG GmbH

X-Ray Protection and Medical
Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДИКА“
Уповноважений
представник в Україні
01014, Україна, м. Київ, вул.
Звіринецька, 63
+38 (044) 228 84 17
+38 (097) 796 71 30
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua