



<b>Models:</b>	<b>RP2761E</b>	Testicle Protection for men
	<b>RP2761K</b>	Testicle Protection for boys
	<b>RP2761G6232</b>	Testicle Shield Set for Children aged 1 to 10

Made in Germany

CE 0302

## RP276 – Testicle Shield for Men/Boys/Children

This male gonad protection model covers the penis and scrotum all around with a lead equivalent of 1.00 mm Pb. Of course, MAVIG has made sure that there is no risk of injury to the patient when using the protection.

The hook-and-loop fastener provides a good hold when adjusted correctly. For cleaning and disinfection, this model series can be unfolded using the pressstuds. We also have suitable disposable hygiene bags\* in our range.

- **Protective Material** Standard lead
- **Lead Equivalent** 1.00 mm Pb
- **Outer Cover/Colour** Soft PVC cover in light grey (boys/men);  
Textile cover made of special, easy to clean material  
in design Happy (infants/baby/children)
- **Sizes**
  - RP2761E - for men (dimensions approx. 155 x 120 x 90 mm)
  - RP2761K - for boys (dimensions approx. 110 x 90 x 65 mm)
  - RP2761G6232 - Set with 3 protections from infant to child sizes

**Hint:** Different hygienic cover sets are available for the sizes E and K, which each contain 100 covers. The set RP-HYG-CEK contains 100 covers for size E and 100 covers for size K.

### MAVIG GmbH

X-Ray Protection and Medical  
Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 796 71 30  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua