



AW417 - Wall Mount

Models:	RP6435S	Gonad apron (toddlers)	RP6435L	Gonad apron (adults)
	RP6435M	Gonad apron (children)	RP6435X	Gonad apron (X-large)

Made in Germany

CE 0302

RP643 – Gonad Apron for Patients

An article for daily use makes special demands on material quality and workmanship but the functionality for user-friendly use must also be allowed for. Even with the high patient flow in today's radiology practices, the RP643 is ideally suited, as the X-ray protection is quick and practical to put on and take off. It couldn't be easier than with the flexible clamping mechanism of our RP643 series. The radiation protection for your patients is perfect and the fit is just right.

In producing our gonad aprons, we use only highquality, flexible standard lead material. The best radiation protection properties are, of course, a basic requirement.

- **Protective Material** Standard lead
- **Lead Equivalent** 0.50 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:
Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy
- **Sizes**
 - S - for small children, covers 20 x 25 cm (H x W)
 - M - for children, covers 30 x 30 cm (H x W)
 - L - for adults, covers 37 x 40 cm (H x W)
 - XL - for adults, covers 45 x 50 cm (H x W)

MAVIG GmbH

X-Ray Protection and Medical Suspension Systems
Stahlgruberring 5
81829 Munich
Germany
PO Box 82 03 62
81803 Munich
Germany
Phone +49 (0) 89 / 420 96 - 0
Fax: + 49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com
www.mavig.com

ТОВ „ОНКОМЕДІКА“ Уповноважений представник в Україні

01014, Україна, м. Київ, вул.
Звіринецька, 63
+38 (044) 228 84 17
+38 (097) 796 71 30
+38 (097) 959 10 89
info@oncomedica.ua
www.oncomedica.ua
www.mavig.com.ua