



EN

RP668

NEW IN DESIGN "HAPPY"

MAVIG



<b>Models:</b>	<b>RP66852S</b>	Children's X-Ray Protective Coat (SMALL)
	<b>RP66852M</b>	Children's X-Ray Protective Coat (MEDIUM)
	<b>RP66852L</b>	Children's X-Ray Protective Coat (LARGE)

Made in Germany

CE 0302

## RP668 – Children's All-round Protective Coat

The cut combines ideal radiation protection with variable fit and high wearing comfort. The flexible material and generous width adjustment allow for a comfortable fit.

Children's X-ray protection, in particular, must be easy to put on and take off. So, you can adjust the RP668 children's coat quickly and securely using the hook-and-loop fasteners at the front. Children's X-ray protection is manufactured by MAVIG with the innovative NovaLite radiation protection material to achieve the best possible weight ratio and thus impose only the most necessary physical strain on the children.

- **Protective Material** Novalite - radiation protective material
- **Lead Equivalent** Front protection 0.50 mm Pb / Back panels 0.25 mm Pb
- **Outer Cover/Colour** ComforTex® HPMF, ComforTex® HPMF Plus (design Happy only), available in the colours:  
*Titan, Blue Admiral, Regatta, Curacao, Lagoon, Mandarin, Orchidee, Baccara, or design: Happy*
- **Size**

<b>SMALL</b>	Length: 60 cm, Hip girth: 68 - 75 cm for children 3 to 5 years old
<b>MEDIUM</b>	Length: 80 cm, Hip girth: 76 - 83 cm for children 6 to 8 years old
<b>LARGE</b>	Length: 84 cm, Hip girth: 84 - 95 cm for children 9 to 12 years old

**MAVIG GmbH**  
X-Ray Protection and Medical  
Suspension Systems  
Stahlgruberring 5  
81829 Munich  
Germany  
PO Box 82 03 62  
81803 Munich  
Germany  
Phone +49 (0) 89 / 420 96 - 0  
Fax: + 49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com  
www.mavig.com  
**ТОВ „ОНКОМЕДІКА“**  
**Уповноважений**  
**представник в Україні**  
01014, Україна, м. Київ, вул.  
Звіринецька, 63  
+38 (044) 228 84 17  
+38 (097) 796 71 30  
+38 (097) 959 10 89  
info@oncomedica.ua  
www.oncomedica.ua  
www.mavig.com.ua