

EN

RP679

MAVIG



**Sets:** **RP6791E12** Ovarian protection with stretch belt RP679 for adults (set with 3 sizes)  
**RP6791K12** Ovarian protection with stretch belt RP679 for children (set with 3 sizes)

Made in Germany

CE 0302

## RP679 – Ovarian Protection with Stretch Belt

The patient shield RP679 protects the ovaries during X-rays of standing patients.

It is also ideally suited for children as patient protection in the appropriate size set for paediatrics and for X-rays in medical diagnostics. The stretch belt\* with snap buckle closure makes it easy to put on and off and ensures convenient size adjustment.

- **X-ray protective material Standard lead with a lead equivalent of 1.00 mm Pb**
- **Comfortable size adjustment by means of the stretch belt with hygienic cover and snap buckle**
- **Quick and easy to put on**
- **Set consisting of the stretch belt\* and three protectors with the following dimensions:**  
 for adults: Small (Blue Admiral, 16 x 11 cm), Medium (Mandarin, 18 x 13 cm), Large (Regatta, 20 x 19 cm)  
 for children: Small (Orchidee, 7 x 6 cm), Medium (Lagoon, 9 x 8 cm), Large (Curacao, 10 x 10 cm)
- **Protectors can be attached to the belt with the hook & burr fastening**

\* A custom size belt is available for patients with a girth greater than 128 cm.

### MAVIG GmbH

X-Ray Protection and  
 Medical Suspension  
 Systems Stahlgruberring 5  
 81829 Munich  
 Germany  
 PO Box 82 03 62  
 81803 Munich  
 Germany  
 Phone +49 (0) 89 / 420 96  
 - 0 Fax: + 49 (0) 89 / 420  
 96 - 200 e-Mail  
 info@mavig.com  
 www.mavig.com

**ТОВ „ОНКОМЕДИКА“**  
**Уповноважений**  
**представник в Україні**  
 01014, Україна, м. Київ,  
 вул. Звіринецька, 63  
 +38 (044) 228 84 17  
 +38 (097) 796 71 30  
 +38 (097) 959 10 89  
 info@oncomedica.ua  
 www.oncomedica.ua  
 www.mavig.com.ua