

Lower Body Protection - UT69

MAVIG



- UT6901** Lower Body Protection with two removable upper shield (25.0 and 17.0 cm)
- UT6902** Lower Body Protection with one removable upper shield (25.0 cm)
- UT6904** Lower Body Protection with one removable upper shield (17.0 cm)



An optional accessories rail allows for an e.g. additional protective panel or for small devices to be hung, if compatible.

CE

Made in Germany

Lower Body Protection – UT69

Lower body protection systems are an integral part of every radiation protection concept for interventional radiology. These systems absorb secondary radiation underneath the table and shields the lower extremities and genital area of the examiner.

The UT69 offers maximum protection, thanks to the intermediate double articulated joint and additional lateral lead panel (removable). The upper shield can be removed quickly in case of emergency or for comfortable patient positioning.

- **Overlapping, flexible X-ray protective panels with PVC covering (light grey)**
- **Lead equivalent Pb 0.50 mm**
- **Suitable for flat profile rails from 7 x 25 mm to 10 x 30 mm**
- **Height: 90.0 cm**
- **Widths: 64.5 cm (main part, 4 overlapping panels)
17.0 cm (additional side panel)**
- **Total width of all lower body protective panels: max. 80.0 cm (overlapping, angle)**
- **Two separate removable upper shields with 25.0 and 17.0 cm height, 60.0 cm width**
- **Wall mounts included for storage**

MAVIG GmbH

X-Ray Protection & System Solutions

Stahlgruberring 5 81829
Munich Germany

Phone: 089 / 420 96 – 0

Fax: 089 / 420 96 – 200

e-Mail: info@mavig.com
www.mavig.com

ТОВ "ОНКОМЕДИКА"

Уповноважений представник в Україні

01014, Україна, м. Київ, вул.

Звіринська, 63

+38 (044) 228 84 17

+38 (097) 959 10 89

info@oncomedica.ua

www.oncomedica.ua

www.mavig.com.ua